



Pressure Redistribution Pads

Reducing Pressure
Preventing Injury
Reducing Costs

dermisplus[®]
PREVENT

Now
Available on
Drug Tariff
& NHS SC

Reducing Pressure. Reducing Costs.

Preventing and managing pressure damage can be complicated, time-consuming and a drain on clinical resources. The use of medical devices can also be a source of potential pressure damage. It is estimated that it costs the NHS £3.8m per day to treat the 700,000 people developing a pressure ulcer each year (NHS Improvement 2016).

To help care providers in the prevention and management of Category 1 pressure damage, Dermisplus Prevent is a new range of pressure redistribution pads designed to help reduce

the risk of pressure damage as part of a pressure area care strategy. They are proven to produce lower peak pressures than competitor product and can provide cost savings of up to 22%.


Available in a range of shapes, sizes and thicknesses, Dermisplus Prevent products are durable and conformable making them ideal for use over bony prominences or with medical devices.

How does it work?

Dermisplus Prevent redistributes peak pressures over a wider surface area so reducing the risk of pressure related tissue damage.

It can be used anywhere on the body which is identified as being 'at risk', especially where there are bony prominences for example, on heels, elbows and sacral areas but also under medical devices such as, medical tubing, mask straps and other equipment.

Pressure redistribution through Dermisplus Prevent



Reducing Pressure

Dermisplus Prevent Offers a 10% Reduction in Peak Pressures¹

The results of comparative testing of Dermisplus Prevent sheets in comparison to Aderma dermal pads, found that Dermisplus Prevent offers on average a **10% better peak pressure redistribution** – so providing more effective pressure redistribution and helping to protect your patients from pressure damage.

Reducing Costs


Using Dermisplus Prevent could reduce spend on pressure redistribution pads by up to 22% versus competitor products and so helping to reduce the financial burden associated with pressure ulcer prevention and management.

How to retain

Dermisplus Prevent products need to be secured in place using an appropriate retention product such as:

- Tape
- Tubular bandage
- Underwear

Peak Pressure Analysis



Indications

Dermisplus Prevent is indicated for use on:

- Intact skin where there is a risk of pressure damage, including from medical devices
- Areas of the body with existing Category 1 pressure damage
- Recently healed ulcers to protect against the risk of re-ulceration

Dermisplus Prevent has **10% better** peak pressure redistribution than competitor products¹

Save up to 22%¹

Product and Application Guide



Benefits



Effective

Dermisplus Prevent reduces peak pressures by 10% in comparison to competitor products



Cost effective

In a comparison of usage for one acute Trust, the use of Dermisplus Prevent could save up to 22%



Durable

Dermisplus Prevent dressings remain stable in use with residue free removal, making dressing changes easier



Flexible Application

Each of the Dermisplus Prevent dressings can be used across many different applications making the range versatile and able to be adapted to individual patient needs. For example, the heel pads can be used on elbows and the sacral dressing around ankles



Easy Application

Dermisplus Prevent dressings are flexible and conformable making them easy to apply



Cut to size

Every patient is different and so Dermisplus Prevent dressings can easily be cut to size without any fibre shedding



Washable

The dressings are designed to be washed and re-used on the same patients. Simply wash with soap and water and leave to dry



Patient friendly

Dermisplus Prevent is latex-free, non-irritating and skin friendly, reducing the risk of allergies and increasing patient comfort

Dermisplus Prevent Product Range Order Codes

Product Code	UOI	Product Description	NPC	PIP Code
8500007	1 x 5	Dermisplus Prevent Pad 10cm x 10cm x 0.3cm	ELY718	4102471
8500003	1 x 5	Dermisplus Prevent Pad 10cm x 10cm x 1.2cm	ELY731	4102489
8500005	1 x 2	Dermisplus Prevent Pad 20cm x 20cm x 0.3cm	ELY720	4102497
8500009	1 x 2	Dermisplus Prevent Pad 20cm x 20cm x 1.2cm	ELY732	4102505
8500006	1 x 5	Dermisplus Prevent Strip 30cm x 5cm x 0.3cm	ELY722	4102513
8500001	1 x 5	Dermisplus Prevent Strip 50cm x 2.5cm x 0.3cm	ELY723	4102521
8500002	1 x 2	Dermisplus Prevent Heel Standard	ELY729	4102547
8500004	1 x 2	Dermisplus Prevent Heel Large	ELY730	4102554
8500008	1 x 1	Dermisplus Prevent Sacrum	ELY724	4102539
8500010	5 x 10s	Dermisplus Prevent Tubes 0.2cm x 10cm	ELY983	4198909

1. Data on file

For more information:

Call +44 (0)330 460 6030

Email: info@frontier-group.co.uk

Website: www.frontier-group.co.uk



Frontier Medical Group, Innova One, Tredegar Business Park, Tredegar, Blaenau Gwent NP22 3EL United Kingdom.

© Copyright Frontier Therapeutics Limited, 2022
All rights reserved. GP0132_A42_09

dermisplus[®]
PREVENT